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$$f = \frac{S \cdot d \cdot b^2}{2 m_1 \cdot t^2 \cdot L} \quad (II)$$

The device can be used to demonstrate the gravitational attractive force between masses and to determine the gravitational constant  $f$ .

**Measurement principle** (Fig. 1)

The static equilibrium (position I) of a torsion pendulum with a period of approx. 10 min. is disturbed by a change in position of the outer masses  $m_1$ , which affect the dumbbell-shaped pendulum body. The oscillations become damped and the pendulum takes up a new equilibrium position (position 2). The angle between the two positions is a measure of the active gravitational force.

The oscillation of the pendulum, which is equipped with a concave mirror, is indicated by a light pointer. This is possible using

- visible light directly on a mm-scale or
- infrared light along with the IR position detector (332 11). The latter enables the measurement values to be plotted or evaluated by computer.

The gravitational constant  $f$  can be obtained from the oscillation curve with respect to time, the mass  $m_1$  and the geometry of the arrangement using either the end deflection method or (in a quicker process) the acceleration method.

In the *end deflection method* the torsion pendulum period  $T$  and the distance  $S$  separating the light pointer positions are evaluated for the two equilibrium positions:

$$f = \frac{\pi^2 \cdot b^2 \cdot d \cdot S}{m_1 \cdot T^2 \cdot L} \quad (I)$$

In the *acceleration method*, the acceleration of the torsion pendulum  $a = 2s/t^2$  is evaluated using the masses  $m_1$  after the equilibrium position has been disturbed:

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## 1 Safety notes

Protect the sensitive bronze band of the torsion pendulum from uncontrolled mechanical loading.

- Locking screws (7) of the oscillating system (see Fig. 2) before the device has been correctly assembled and brought into position.
- Always lock the oscillating system when the device is not in use; in particular, make sure it is locked during transport and assembly.

The knurled screw (3.1) - see Fig. 2 - for fixing the torsion head (pre-adjusted on delivery) is only to be slightly loosened if the adjustment of the zero point proves necessary on putting the correctly assembled torsion balance into operation. The grub screw (3.2) - see Fig. 2 - for fixing the pendulum holder is only to be loosened when replacing the torsion band as described in section 4.

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## 2 Description, technical data, scope of supply (see Fig. 2)


Fig.2  
Gravitations-Drehwaage.  
Gravitation torsion balance

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- ① Metal housing (15 cm dia.), with metal ring and glass covers; two-part slide (1.1) to prevent interference from convection currents between the housing and protective tube ②
- ② Protective tube (25 cm long) for torsion band ④
- ③ Torsion head with pendulum holder, can be rotated for adjustment of the pendulum equilibrium position when screw has been loosened (3.1); pendulum holder fixed with grub screw (3.2)  
**Important!**  
**Do not loosen the knurled screw (3.1) unless it is essential to carry out a zero-point adjustment in accordance with section 3.5. Only loosen the grub screw (3.2) when the torsion band is being replaced as described in section 4.**
- ④ Torsion band made from bronze, 26 cm long  
Spare part no.: 683 21
- ⑤ Dumbbell-shaped pendulum body, consisting of 2 lead balls (5.1) on a metal rod

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- ⑥ Concave mirror for light pointer indication of pendulum motion  
Focal length  $f$  approx. 30 cm
- ⑦ Screws for locking the pendulum body ⑤ using pair of springs (7.1)
- ⑧ Pairs of large lead balls  
Spare part no. for 1 ball: 683 22
- ⑨ Ball carrier, can be rotated around stand rod ⑪, for correct experimental arrangement of the large lead balls ⑧
- ⑩ Supporting ring with fixing screw (10.1) for ball carrier
- ⑪ Stand rod (9 cm x 1.2 cm dia.) for assembly of the device in stand material

Included in scope of supply:  
1 m self-adhesive scale tape with cm and mm divisions

### Important experimental data (see Fig. 1):

Housing depth: 30 mm

Torsion pendulum

Period: approx. 10 min

Diameter / mass  $m_2$  of a lead ball (5.1):  
15 mm / 20 g

Distance  $d$  between the center of a ball and the axis of rotation: 50 mm

Diameter / mass  $m_1$  of a large lead ball ⑥:  
64 mm / 1.5 kg  $\pm 5$  g

Distance  $b$  between the center of the large ball (when contact is made with housing) and the small ball (in the equilibrium position): 47 mm

## 3 Operation

### Important!

Satisfactory experiment results are only possible when the torsion pendulum has been properly adjusted and the torsion oscillations produced by attraction between the masses are not affected by unwanted pendulum movements. The pendulum is very sensitive to any disturbance of the experiment setup: make sure that the experiment setup is absolutely stable, e.g. by attaching it to a solid wall or placing it on a sturdy bench or table.

Temperature variations cause convection in the housing of the torsion balance, which in turn cause undesired motions of the torsion pendulum (see Fig. 6).

For this reason, select an experiment site which does not stand in direct sunlight or drafts.

### 3.1 Additional equipment required

#### 3.1.1. Oscillation indication using (visible) light marks on a mm-scale

1 Lamp housing	450 60
1 Lamp, 6 V, 30 W	450 51
1 Spherical condenser (with 1-mm slit screen)	460 20
1 A.C. voltage source, 6 V, 30 W e.g.	562 73

✓ Steel tape measure	311 77
✓ Stop-clock e.g.	313 05

1 Suspended magnet	510 44
--------------------	--------

#### Wall assembly (see Fig. 3)

1 Large stand base	300 31
1 Pair of levelling screws	300 06
1 Leybold multiclamp	301 01
1 Rotatable clamp	301 03
1 Stand rod, 47 cm	300 42

Assembly aids and accessories:

Electric drill, masonry drill bit, wall plugs (6 mm)

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Setup for the gravitation torsion balance; display of the oscillation curve by means of (visible) light pointer on mm-scale (not suitable for IR position detector 332 11)

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Benchtop assembly of the gravitation torsion balance; electronic determination of the oscillation curve using an IR position detector (332 11); recording with TY-recorder or further processing with computer assistance (see Fig. 4.1)

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3.1.2 Recorder plotting or computer-assisted recording of the oscillation using infrared light; rail assembly on a bench (see Fig. 4)

1 Infrarot-Position-Detector 332 11  
1 Wechselspannungsquelle, 12 V z.B. 562 73

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1 TY-recorder e.g. 575 701  
or  
MS DOS computer via serial interface  
RS 232 cable e.g. 530 008  
Program with experiment-specific default settings e.g. contained on the demo disk for "Universal Data Acquisition", included in the scope of supply of the IR position detector (332 11); see measuring example Fig. 4.1.

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or  
MS DOS computer with interface  
CASSYpack-E 524 007  
Program "Measuring and Evaluating" 524 112

1 Rollbandmaß 311 77  
1 Schwebemagnet 510 44

1 Steel tape measure 311 77  
1 Suspended magnet 510 44

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**Note:**

Further information on the use of the infrared position detector can be found in the corresponding instruction sheet.

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Computer-assisted recording of the oscillations of the gravitation torsion balance

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### 3.2 Mounting the stand base on the wall (required for wall mounting as shown in Fig. 3 before using the setup for the first time.)

#### Important:

Wall mounting of the stand material is not suitable for the experiment setup with the IR position detector 332 11 (see section 3.1.2; Fig. 4).

Correctly mount the stand base on its apex hole as shown in Fig. 5.1 - 3.3 and align parallel to the wall using levelling screws; use solid base material if the levelling screws press into the wall.


Fig. 5.1-5.3

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### 3.3 Preparing for the experiment

Assemble the arrangement - initially without the large lead balls - on the wall as shown in Fig. 3 (or similarly on the optical bench on a sturdy lab bench).

Align the lamp housing so that the reflected beam is projected unobstructed onto the scale, which should be positioned at least 5 m distant.

Align the filament so that it is vertical before attaching the condenser with screen holder and slit (vertically!) to the lamp housing.

Move the insert in the lamp housing so that a focused image of the filament is projected onto the torsion balance mirror. To do this, hold a white sheet of paper, on which the filament forms an image, directly in front of the balance.

Move the lamp housing so that the slit is focused onto the scale. If necessary, adjust the zero-point as described in section 3.5.

Attach the lead balls and bring the arrangement to an extreme position. In carrying out this step make absolutely sure that the housing does not come into contact with either your finger or a lead ball.

### 3.4 Experiment procedure; measuring example

#### Important!

After the experiment has been set up as described in section 3.3, it must be left undisturbed for at least two hours so that the pendulum can settle in its equilibrium state.

Turn the ball carrier without touching the housing.

To shorten this settling time, make use of the diamagnetism of the lead: each time one of the pendulum lead balls approaches the glass cover, hold up a strong magnet - without touching the glass - until the motion of the system is reversed.

The experiment setup must not be subjected to any shocks or temperature variations which may cause convection in the torsion balance housing while recording measured values (see Fig. 6)

Check the zero-point stability before starting the measurement. If necessary, observe and note the zero-point fluctuations for at least 10 minutes so that you can determine an average value for  $x_0$ .

At time  $t = 0$  rotate the carrier with the lead balls quickly from one extreme position to the other. However, in carrying out this movement, make sure that fingers or lead balls do not come into contact with the housing. Start the stop-clock immediately after shifting the carrier.

Read off and note the position of the light pointer on the scale every 30 s (Fig. 6) for at least 3 periods when using the end displacement method and for 1 period when using the acceleration method.

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24-hour log of the position of the torsion pendulum: position changes are essentially due to (day)time dependent variations in temperature (particularly noticeable for exposure to direct sunlight)

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Oscillations of the gravitation torsion balance around the final equilibrium position  $x_{00}$

#### Measuring example:

Equipment constants:

Mass of the large lead ball:  $m_1 = 1.5 \text{ kg}$   
 Distance between ball center and axis of rotation:  $d = 0.05 \text{ m}$   
 Distance between center point of the large ball (when touching housing) and the small ball (in equilibrium position):  $b = 0.047 \text{ mm}$

Distance between mirror and torsion balance:  $L = 4.425 \text{ m}$

Oscillation period  $T$  (from Fig. 7):

$$T = (2790 - 315)/4 \text{ s} = 618.8 \text{ s}$$

Equilibrium position at start:  $x_0 = 47 \text{ cm}$

Equilibrium position at end  $x_{00}$  (determined from three sequential extremes) e.g.

$$x_{00} = \frac{(x_1 + x_3)/2 + x_2}{2} = \frac{x_1}{4} + \frac{x_2}{2} + \frac{x_3}{4} = 62.3 \text{ cm}$$

Difference  $S$  between light pointer positions for the initial and final pendulum equilibrium states:

$$S = x_{00} - x_0 = 62.3 \text{ cm} - 47 \text{ cm} = 15.3 \text{ cm}$$

Substitution in equation (1) for the end displacement method gives the gravitation constant  $f$  as

$$f = \frac{\pi^2 \cdot (0.047 \text{ m})^2 \cdot 0.05 \text{ m} \cdot 0.153 \text{ m}}{1.5 \text{ kg} \cdot (618.8 \text{ s})^2 \cdot 4.425 \text{ m}} = 6.56 \cdot 10^{-11} \text{ m}^3 \text{ kg}^{-1} \text{ s}^2$$

### 3.5 Zero-point adjustment

The objective of zero-point adjustment is to position the torsion pendulum so that the dumbbell-shaped pendulum body ⑤ (see Fig. 2) is aligned parallel to the glass plates of the housing in the equilibrium position.

An adjustment is necessary:

possibly before first using the system if the manufacturer's adjustment has been modified, e.g. as a result of transportation, after improper handling (through uncontrolled rotation of torsion head) ③,

after a new torsion band has been inserted as described in section 4.

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#### **Preparation for zero-point adjustment:**

Assemble the torsion balance with light pointer arrangement on the lab bench e.g. as described in section 3.3 (or in accordance with the corresponding information in the Instruction Sheet 332 11), but without the large lead balls.

Unlock the torsion pendulum by loosening the screws ⑦ (see Fig. 2).

Using the levelling screws on the stand base (or on the optical bench), align the torsion balance so that the torsion pendulum hangs vertically (check: pin at the end of the pendulum hangs in the axis of the hole).

After these preparatory steps have been carried out, allow the arrangement to hang disturbance-free for around 1 day and then carry out the zero-point adjustment.

#### **Zero-point adjustment:**

Lock the torsion pendulum with screw ⑦ and set up the torsion balance at the experiment location (e.g. wall mounting as shown in Fig. 3); then unlock the pendulum to set it in motion. Note the end deflections of the light pointer on both sides of the scale.

Observe the equilibrium position to which the system tends; to shorten the observation time damp the oscillations diamagnetically: each time one of the pendulum lead balls approaches the glass cover, use a strong magnet to reverse the motion of the system - without touching the glass.

Should the light pointer not show a tendency towards an equilibrium position roughly in the center between the two noted direction reversal points, unscrew the knurled screws (3.1) and turn the torsion head ③ through a small angle in the direction of the "desired zero-point".

Once again observe the two end deflections and the zero-point tendency of the light pointer.

Continue this zero-point adjustment until the "desired zero-point" in the middle between the two direction reversal points of the light pointer has been reached.

## **4 Replacing the torsion band**

A defective torsion band can be replaced by a new band which is supplied in a ready-for-assembly state as spare part (683 21). The replacement is carried out as follows:

#### **Preparation**

- Assemble the torsion balance (with locked pendulum body) on a stand base which can be levelled (e.g. 300 01/06) in such a way that the back side is at the front; to ensure safe handling, support the protective tube with replacement torsion band (683 21) in a clamp (Fig. 7).
- Remove screw (a<sub>1</sub>); subsequently loosen the other screws fixing the glass plate until it can be removed (Fig. 8.2).
- Remove glass plate (a<sub>2</sub>) and cardboard seal (a<sub>3</sub>) and pull out the front section of the cardboard slide (a<sub>4</sub>).

#### **Removing the defective torsion band**

- Loosen grub screw (b<sub>1</sub>) while holding the band holder (b) tightly. Remove the band holder (Fig. 8.3).
- If the torsion band is broken, remove end piece (c) (Fig. 8.4).

#### **Inserting the new torsion band**

- Loosen screw x<sub>1</sub> and remove screw x<sub>2</sub> while holding band holder (B) tightly.
- Very carefully remove the arrangement from the protective cover and "thread" it into the torsion balance until the end piece (C) hangs 1 mm to 2 mm above the hole (c<sub>2</sub>). Make sure the very sensitive torsion band does not come into contact with the walls. Subsequently fix pendulum holder (B) by tightening screw ((b<sub>1</sub>) (Fig. 8.6).
- Wait until the torsion band has settled and, if necessary, use the levelling screws on the stand base to align the torsion band so that the end piece (C) hangs vertically above the hole.

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- Loosen screw (b<sub>1</sub>). Hold the band holder (B) tightly and lower it a few mm - without twisting it - until the machined part of the end piece (C) has been completely guided into the hole (c<sub>1</sub>); do not release the band holder until it has been carefully fixed in place by tightening the grub screw (b<sub>1</sub>) (Fig. 8.7); finally fix the end piece (C) with screw (c<sub>1</sub>).
- Insert the front section of the slide (a<sub>4</sub>). After correctly reinserting the cardboard disc (a<sub>3</sub>) and the glass plate (a<sub>2</sub>), close the housing with screw (a<sub>1</sub>).
- Insert cardboard slide (a<sub>4</sub>), and then reclose the housing with screw (a<sub>1</sub>) after correctly inserting sealing disc (a<sub>3</sub>) and glass plate (a<sub>2</sub>).

After assembly has been carried out, unlock the torsion pendulum and check if the pendulum body can oscillate freely; should it lie on the locking springs despite unlocking, slightly raise the band holder (B), with screw (b<sub>1</sub>) loosened, and again carefully fix in place.

The torsion pendulum must hang undisturbed for at least 12 hours before the required zero-point adjustment as described in section 3.5 can be carried out.


Fig. 8.1


Fig. 8.3


Fig. 8.5


Fig. 8.2


Fig. 8.4


Fig. 8.7


Fig. 8.6

The enormous strength of the Earth's attraction for the small masses, in comparison with their attraction for the large masses, is what originally made the measurement of the gravitational constant such a difficult task. The torsion balance (invented by Charles Coulomb) provides a means of negating the otherwise overwhelming effects of the Earth's attraction in this experiment. It also provides a force delicate enough to counterbalance the tiny gravitational force that exists between the large and small masses. This force is provided by twisting a very thin beryllium copper ribbon.

The large masses are first arranged in Position I, as shown in Figure 12, and the balance is allowed to come to equilibrium. The swivel support that holds the large masses is then rotated, so the large masses are moved to Position II, forcing the system into disequilibrium. The resulting oscillatory rotation of the system is then observed by watching the movement of the light spot on the scale, as the light beam is deflected by the mirror.

Any of three methods can be used to determine the gravitational constant,  $G$ , from the motion of the small masses. In Method I, the final deflection method, the motion is allowed to come to resting equilibrium—a process that requires several hours—and the result is accurate to within approximately 5%. In method II, the equilibrium method, the experiment takes 90 minutes or more and produces an accuracy of approximately 5% when graphical analysis is used in the procedure. In Method III, the acceleration method, the motion is observed for only 5 minutes, and the result is accurate to within approximately 15%.

### METHOD I: Measurement by Final Deflection

Setup Time: ~ 45 minutes; Experiment Time: several hours  
Accuracy: ~ 5%

#### Theory

With the large masses in Position I (Figure 13), the gravitational attraction,  $F$ , between each small mass ( $m_2$ ) and its neighboring large mass ( $m_1$ ) is given by the law of universal gravitation:

$$F = Gm_1m_2/b^2 \quad (1.1)$$

where  $b$  = the distance between the centers of the two masses.

► Note: 5% accuracy is possible in Method I if the experiment is set up on a sturdy table in an isolated location where it will not be disturbed by vibration or air movement.

► Note: 5% accuracy is possible in Method II if the resting equilibrium points are determined using a graphical analysis program.


Figure 13  
Origin of variables  $b$  and  $d$

The gravitational attraction between the two small masses and their neighboring large masses produces a net torque ( $\tau_{grav}$ ) on the system:

$$\tau_{grav} = 2Fd \quad (1.2)$$

where  $d$  is the length of the lever arm of the pendulum bob crosspiece.

Since the system is in equilibrium, the twisted torsion band must be supplying an equal and opposite torque. This torque ( $\tau_{band}$ ) is equal to the torsion constant for the band ( $\kappa$ ) times the angle through which it is twisted ( $\theta$ ), or:

$$\tau_{band} = -\kappa\theta. \quad (1.3)$$

Combining equations 1.1, 1.2, and 1.3, and taking into account that  $\tau_{grav} = -\tau_{band}$ , gives:

$$\kappa\theta = 2dGm_1m_2/b^2$$

Rearranging this equation gives an expression for  $G$ :

$$G = \frac{\kappa\theta b^2}{2dm_1m_2} \quad (1.4)$$

To determine the values of  $\theta$  and  $\kappa$  — the only unknowns in equation 1.4 — it is necessary to observe the oscillations of the small mass system when the equilibrium is disturbed. To disturb the equilibrium (from  $S_1$ ), the swivel support is rotated so the large masses are moved to Position II. The system will then oscillate until it finally slows down and comes to rest at a new equilibrium position ( $S_2$ ) (Figure 14).

At the new equilibrium position  $S_2$ , the torsion wire will still be twisted through an angle  $\theta$ , but in the opposite direction of its twist in Position I, so the total change in angle is equal to  $2\theta$ . Taking into account that the angle is also doubled upon reflection from the mirror (Figure 15):

$$\begin{aligned} \Delta S &= S_2 - S_1, \\ 4\theta &= \Delta S/L \quad \text{or} \\ \theta &= \Delta S/4L \end{aligned} \quad (1.5)$$


Figure 14  
Graph of Small Mass Oscillations


Figure 15  
Diagram of the experiment showing the optical lever.

The torsion constant can be determined by observing the period ( $T$ ) of the oscillations, and then using the equation:

$$T^2 = 4\pi^2 I / \kappa \quad (1.6)$$

where  $I$  is the moment of inertia of the small mass system.

The moment of inertia for the mirror and support system for the small masses is negligibly small compared to that of the masses themselves, so the total inertia can be expressed as:

$$I = 2m_2(d^2 + 2/5r^2) \quad (1.7)$$

Therefore:

$$\kappa = 8\pi^2 m_2 \frac{d^2 + 2/5r^2}{T^2} \quad (1.8)$$

Substituting equations 1.5 and 1.8 into equation 1.4 gives:

$$G = \pi^2 \Delta S b^2 \frac{(d^2 + 2/5r^2)}{T^2 m_1 L d} \quad (1.9)$$

All the variables on the right side of equation 1.9 are known or measurable:

$$r = 9.55 \text{ mm}$$

$$d = 50 \text{ mm}$$

$$b = 46.5 \text{ mm}$$

$$m_1 = 1.5 \text{ kg}$$

$$L = (\text{Measure as in step 1 of the setup.})$$

By measuring the total deflection of the light spot ( $\Delta S$ ) and the period of oscillation ( $T$ ), the value of  $G$  can therefore be determined.

## Procedure

1. Once the steps for leveling, aligning, and setup have been completed (with the large masses in Position I), allow the pendulum to stop oscillating.
2. Turn on the laser and observe the Position I end point of the balance for several minutes to be sure the system is at equilibrium. Record the Position I end point ( $S_I$ ) as accurately as possible, and indicate any variation over time as part of your margin of error in the measurement.

- Carefully rotate the swivel support so that the large masses are moved to Position II. The spheres should be just touching the case, but take care to avoid knocking the case and disturbing the system.

**Note:** You can reduce the amount of time the pendulum requires to move to equilibrium by moving the large masses in a two-step process: first move the large masses and support to an intermediate position that is in the midpoint of the total arc (Figure 16), and wait until the light beam has moved as far as it will go in the period; then move the sphere across the second half of the arc until the large mass support just touches the case. Use a slow, smooth motion, and avoid hitting the case when moving the mass support.


- Immediately after rotating the swivel support, observe the light spot and record its position ( $S_1$ ).
- Use a stop watch to determine the time required for one period of oscillation ( $T$ ). For greater accuracy, include several periods, and then find the average time required for one period of oscillation.

**Note:** The accuracy of this period value ( $T$ ) is very important, since the  $T$  is squared in the calculation of  $G$ .

- Wait until the oscillations stop, and record the resting equilibrium point ( $S_2$ ).

### Analysis

- Use your results and equation 1.9 to determine the value of  $G$ .
- The value calculated in step 2 is subject to the following systematic error. The small sphere is attracted not only to its neighboring large sphere, but also to the more distant large sphere, though with a much smaller force. The

**Figure 16**  
Two-step process of moving the large masses to reduce the time required to stop oscillating

geometry for this second force is shown in Figure 17 (the vector arrows shown are not proportional to the actual forces).

From Figure 17,

$$f = F_0 \sin \Phi$$

$$\sin \Phi = \frac{b}{(b^2 + 4d^2)^{1/2}}$$

The force,  $F_0$  is given by the gravitational law, which translates, in this case, to:

$$F_0 = \frac{Gm_2m_1}{(b^2 + 4d^2)}$$

and has a component  $f$  that is opposite to the direction of the force  $F$ :

$$f = \frac{Gm_2m_1b}{(b^2 + 4d^2)(b^2 + 4d^2)^{1/2}} = \beta F$$

This equation defines a dimensionless parameter,  $b$ , that is equal to the ratio of the magnitude of  $f$  to that of  $F$ . Using the equation  $F = Gm_1m_2/b^2$ , it can be determined that:

$$b = b^3/(b^2 + 4d^2)^{3/2}$$

From Figure 17,

$$F_{net} = F - f = F - bF = F(1 - b)$$

where  $F_{net}$  is the value of the force acting on each small sphere from *both* large masses, and  $F$  is the force of attraction to the nearest large mass only.

Similarly,

$$G = G_0(1 - b)$$

where  $G$  is your experimentally determined value for the gravitational constant, and  $G_0$  is corrected to account for the systematic error.

Finally,

$$G_0 = G/(1 - b)$$

Use this equation with equation 1.9 to adjust your measured value.


Figure 17  
Correcting the measured value of  $G$

## METHOD II: Measurement by Equilibrium Positions

Observation Time: ~ 90+ minutes

Accuracy: ~ 5 %

### Theory

When the large masses are placed on the swivel support and moved to either Position I or Position II, the torsion balance oscillates for a time before coming to rest at a new equilibrium position. This oscillation can be described by a damped sine wave with an offset, where the value of the offset represents the equilibrium point for the balance. By finding the equilibrium point for both Position I and Position II and taking the difference, the value of  $\Delta S$  can be obtained. The remainder of the theory is identical to that described in Method I.

### Procedure


1. Set up the experiment following steps 1-3 of Method I.
2. Immediately after rotating the swivel support to Position II, observe the light spot. Record the position of the light spot ( $S$ ) and the time ( $t$ ) every 15 seconds. Continue recording the position and time for about 45 minutes.
3. Rotate the swivel support to Position I. Repeat the procedure described in step 2.

**Note:** Although it is not imperative that step 3 be performed immediately after step 2, it is a good idea to proceed with it as soon as possible in order to minimize the risk that the system will be disturbed between the two measurements. Waiting more than a day to perform step 3 is not advised.

### Analysis

1. Construct a graph of light spot position versus time for both Position I and Position II. You will now have a graph similar to Figure 18.
2. Find the equilibrium point for each configuration by analyzing the corresponding graphs using graphical analysis to extrapolate the resting equilibrium points  $S_1$  and  $S_2$  (the equilibrium point will be the center line about which the oscillation occurs). Find the difference between

**Note:** To obtain an accuracy of 5% with this method, it is important to use graphical analysis of the position and time data to extrapolate the resting equilibrium positions,  $S_1$  and  $S_2$ .

**Figure 18**  
Typical pendulum oscillation pattern showing equilibrium positions

the two equilibrium positions and record the result as  $\Delta S$ .

3. Determine the period of the oscillations of the small mass system by analyzing the two graphs. Each graph will produce a slightly different result. Average these results and record the answer as  $T$ .
4. Use your results and equation 1.9 to determine the value of  $G$ .
5. The value calculated in step 4 is subject to the same systematic error as described in Method I. Perform the correction procedure described in that section (*Analysis, step 3*) to find the value of  $G_0$ .

### **METHOD III: Measurement by Acceleration**

Observation Time: ~ 5 minutes

Accuracy: ~ 15%

#### **Theory**

With the large masses in Position I, the gravitational attraction,  $F$ , between each small mass ( $m_2$ ) and its neighboring large mass ( $m_1$ ) is given by the law of universal gravitation:

$$F = Gm_1m_2/b^2 \quad (3.1)$$

This force is balanced by a torque from the twisted torsion ribbon, so that the system is in equilibrium. The angle of twist,  $\theta$ , is measured by noting the position of the light spot where the reflected beam strikes the scale. This position is carefully noted, and then the large masses are moved to Position II. The position change of the large masses disturbs the equilibrium of the system, which will now oscillate until friction slows it down to a new equilibrium position.

Since the period of oscillation of the small masses is long (approximately 10 minutes), they do not move significantly when the large masses are first moved from Position I to Position II. Because of the symmetry of the setup, the large masses exert the same gravitational force on the small masses as they did in Position I, but now in the opposite direction. Since the equilibrating force from the torsion band has not

changed, the total force ( $F_{total}$ ) that is now acting to accelerate the small masses is equal to twice the original gravitational force from the large masses, or:

$$F_{total} = 2F = 2Gm_1m_2/b^2 \quad (3.2)$$

Each small mass is therefore accelerated toward its neighboring large mass, with an initial acceleration ( $a_0$ ) that is expressed in the equation:

$$m_2a_0 = 2Gm_1m_2/b^2 \quad (3.3)$$

Of course, as the small masses begin to move, the torsion ribbon becomes more and more relaxed so that the force decreases and their acceleration is reduced. If the system is observed over a relatively long period of time, as in Method I, it will be seen to oscillate. If, however, the acceleration of the small masses can be measured before the torque from the torsion ribbon changes appreciably, equation 3.3 can be used to determine  $G$ . Given the nature of the motion—damped harmonic—the initial acceleration is constant to within about 5% in the first one tenth of an oscillation. Reasonably good results can therefore be obtained if the acceleration is measured in the first minute after rearranging the large masses, and the following relationship is used:

$$G = b^2a_0/2m_1 \quad (3.4)$$


The acceleration is measured by observing the displacement of the light spot on the screen. If, as is shown in Figure 19:

- $\Delta s$  = the linear displacement of the small masses,
- $d$  = the distance from the center of mass of the small masses to the axis of rotation of the torsion balance,
- $\Delta S$  = the displacement of the light spot on the screen, and
- $L$  = the distance of the scale from the mirror of the balance,

then, taking into account the doubling of the angle on reflection,

$$\Delta S = \Delta s(2L/d) \quad (3.5)$$

Using the equation of motion for an object with a constant

**Figure 19**  
Source of data for calculations in Method III

acceleration ( $x = 1/2 at^2$ ), the acceleration can be calculated:

$$a_0 = 2\Delta s/t^2 = \Delta Sd/t^2L \quad (3.6)$$

By monitoring the motion of the light spot over time, the acceleration can be determined using equation 3.6, and the gravitational constant can then be determined using equation 3.4.

### Procedure

1. Begin the experiment by completing steps 1–3 of the procedure detailed in Method I.
2. Immediately after rotating the swivel support, observe the light spot. Record the position of the light spot ( $S$ ) and the time ( $t$ ) every 15 seconds for about two minutes.

### Analysis

1. Construct a graph of light spot displacement ( $\Delta S = S - S_0$ ) versus time squared ( $t^2$ ), with  $t^2$  on the horizontal axis (Figure 20). Draw a best-fit line through the observed data points over the first minute of observation.
2. Determine the slope of your best-fit line.
3. Use equations 3.4 and 3.6 to determine the gravitational constant.
4. The value calculated in step 3 is subject to a systematic error. The small sphere is attracted not only to its neighboring large sphere, but also to the more distant large sphere, although with a much smaller force. Use the procedure detailed in Method I (*Analysis, step 3*) to correct for this force.


Figure 20  
Sample data and best-fit line